

Alcon® Clinical Assessment of a New Multi-Purpose Disinfecting Solution in Asymptomatic and Symptomatic Patients

John Schachet, OD; Lamar Ziegler, OD; David Wakabayashi, OD; Stephen Cohen, OD

ABSTRACT

Purpose: Clinically assess a new multi-purpose disinfecting solution (MPS) with a variety of lens materials in asymptomatic and symptomatic subjects.

Method: Five randomized, double-masked studies conducted at 66 sites evaluated the clinical performance of a new MPS with traditional and silicone hydrogel (SH) soft contact lenses. Regimen 1 (OPTIFREE® Replenish MPS) was compared with Regimens 2 (ReNu® MultiCare™ MPS), 3 (ReNu® MultiCare™ and 4 (Complete® MoistureLUX™). Corneal staining was evaluated by investigators. Subjects recorded rewetting drop use, as well as subjective comfort and symptoms using visual analog scales and Likert-style questionnaires. At a subset of sites, worn lenses were collected for total residual lens lysozyme and lens wettability analyses.

Results: Residual lysozyme of Group IV lenses was significantly lower for Regimen 1 (564.5 and 677.9 µg) compared with Regimen 2 (1002.8 µg) and 3 (692.2 µg, p<0.0001). There was significantly less severe corneal staining with Regimen 1 compared to Regimen 2 in SH lens wearers (p=0.0019) and compared to Regimen 3 in symptomatic Group IV lens wearers (p=0.0089). Group I and IV soft contact lenses cared for with Regimen 1 (58.0° and 58°) had significantly lower contact angles (better wettability) than those using Regimen 2 (76.2° and 90.5°; p<0.0002). Group IV wettability was also significantly better with Regimen 1 (5.3°) compared to Regimen 3 (91.3°; p<0.0001). Rewetting drop use was significantly higher in SH lens wearers using Regimen 4 compared to Regimen 1 (p<0.004). Many indices of comfort were significantly higher for those subjects using Regimen 1, particularly for symptomatic subjects who reported significantly better comfort on all 12 of the Likert statements (p<0.05).

Conclusions: A new MPS was shown to have clinical benefits with traditional and SH lenses when compared with three MPS. This group of studies shows that the choice of lens care regimen can influence ocular health and patient comfort.

INTRODUCTION

With the new lens materials continuously introduced to the market, developing solutions with properties that are beneficial for all soft lenses is an increasing challenge for manufacturers of contact lens care solutions. Solutions need to be able to effectively clean, condition and disinfect lenses, while providing a comfortable lens wearing experience and without causing any negative interactions with the ocular surface. A new multi-purpose disinfecting solution (OPTIFREE® Replenish MPS, Alcon) that reconditions the lens surface for enhanced comfort has been developed and recently marketed. The solution contains the proprietary reconditioning system, TearGlyde™, a combination of the surface active agent IETRONIC® 1304 (registered trademark of BASF Corporation) and a novel surface-active wetting agent, C9-EB3A (nonanoyl ethylenediaminetetraacetate). These components work synergistically with each patient's own natural tears to keep the lens surface moist throughout the day. The new MPS solution contains the proven disinfecting agents, POLYQUAD® (polyquaternium-1) and ADOX® (methylsulfonyl dimethylamine), which have high antibacterial and antifungal activity. This solution also contains propylene glycol, a demulcent and osmolytic agent intended to enhance lens wettability, lubrication and comfort. An extensive clinical effort was undertaken in order to evaluate the safety and efficacy of the new solution. Studies examined the product with various lens materials, including Group I and IV and silicone hydrogel materials, with both asymptomatic and symptomatic subjects, and compared to several multi-purpose solutions.

The purpose of these studies was to clinically assess the safety and efficacy of a new multi-purpose disinfecting solution with a variety of lens materials in both asymptomatic and symptomatic contact lens wearers compared to several lens care regimens.

OBJECTIVE

METHODS

Table 1
Study Summary

Study	Study #	Subjects	Duration	Contact Lenses	Comparison(s)
Study 1	12	285	3 months (Days 0, 7, 14, 21, 28, 30, 30)	Regimen IV	Regimen 1 (OPTIFREE® Replenish MPS)
Study 2	12	233	1 month (Days 0, 7, 14, 30)	Silicone Hydrogel	Regimens 2 & 3 (ReNu® MultiCare™ and Complete® MoistureLUX™)
Study 3	19	343	1 month (Days 0, 7, 14, 28)	Regimen IV	Regimen 3 (Complete® MoistureLUX™)
Study 4	1	80	1 month (Days 0, 1, 4, 30)	Regimen IV	Regimen 2 (ReNu® MultiCare™)
Study 5	22	317	1 month (Days 0, 1, 4, 28)	Regimen IV	Regimen 3 (Complete® MoistureLUX™)

Table 2
Formulation Comparison

Regimen	Cleansers	Antibacterials	Bifens	Additional Wetting Agents
Regimen 1 (OPTIFREE® Replenish MPS)	Sodium Dodecyl Sulfate, C9-EB3A	POLYQUAD®	Sodium Citrate, Sodium Chloride	Propylene Glycol, Glycerol, Polyethylene Glycol, Tetrahydrofuran, Sorbitol, Sorbitan Monostearate, PEG-400, PEG-4000
Regimen 2 (ReNu® MultiCare™)	Propylene Glycol, Polyethylene Glycol, Tetrahydrofuran, Sorbitol, Sorbitan Monostearate, PEG-400, PEG-4000	DIMAC®	Sodium Chloride, Sodium Citrate, Sodium Borate	Propylene Glycol, Glycerol, Polyethylene Glycol, Tetrahydrofuran, Sorbitol, Sorbitan Monostearate, PEG-400, PEG-4000
Regimen 3 (Complete® MoistureLUX™)	Propylene Glycol, Polyethylene Glycol, Tetrahydrofuran, Sorbitol, Sorbitan Monostearate, PEG-400, PEG-4000	Aldehyde	Boic Acid, Sodium Phosphate, Sodium Chloride	Methylcellulose, Polyethylene Glycol, Tetrahydrofuran, Sorbitol, Sorbitan Monostearate, PEG-400, PEG-4000
Regimen 4 (Complete® MoistureLUX™)	Propylene Glycol, Polyethylene Glycol, Tetrahydrofuran, Sorbitol, Sorbitan Monostearate, PEG-400, PEG-4000	Polyquaternium-1, Polyoxyethylene-9, Polyoxyethylene-20, Polyoxyethylene-40, Polyoxyethylene-150	Boic Acid, Sodium Phosphate, Sodium Chloride	Methylcellulose, Polyethylene Glycol, Tetrahydrofuran, Sorbitol, Sorbitan Monostearate, PEG-400, PEG-4000

Table 3
Safety and Efficacy Variables

Variable	Method	Study 1	Unit Group Assessed
Corneal Staining	Visual Questionnaire (1 = severity degree to 5 = severe)	1, 2, 1, 3	All
General Staining (Severity/Area)	Severity (1 = normal, 2 = mild, 3 = moderate, 4 = severe), Area (1 = normal, 2 = mild, 3 = moderate, 4 = severe)	2, 1, 8	All
Wettability	Ex vivo Contact Lens Angle Measurement	4, 5	Regimen 1 (Day 0, 7, 14, 28, 30), Regimen 2 (Day 0, 7, 14, 28, 30)
Lens Cleanliness	Residual Lens Lysozyme (µg) Measurement	1, 5	Day 28 (Study 1) and Day 90 (Study 2)
Rewetting Drop Use	Modified Fluid Overlap Test (data not shown)	All	All
Average Number of Hours per Day	Subjective	1-5	All
Visual Acuity	Investigator-assessed - normal vision, near-normal vision, moderate visual impairment, severe visual impairment, reported by subjects or observed by investigator	1-5	All
Adverse Events	Investigator-assessed	1-5	Post Day 0

RESULTS

Comfort/Acceptability/Symptoms

Figure 1 - Subject Likert Questionnaire
Significant Statistical Differences in Favor of Regimen 1

Statement	Day 14	Day 28
Overall comfort	0.99*	0.81*
Overall acceptability	0.88*	0.77*
Overall eye irritation	0.85*	0.83*
Overall eye redness	0.83*	0.83*
Overall eye itching	0.81*	0.81*
Overall eye watering	0.81*	0.81*
Overall eye burning	0.80*	0.80*
Overall eye dryness	0.79*	0.79*
Overall eye soreness	0.78*	0.78*
Overall eye stinging	0.77*	0.77*
Overall eye irritation	0.76*	0.76*
Overall eye redness	0.75*	0.75*
Overall eye itching	0.74*	0.74*
Overall eye watering	0.73*	0.73*
Overall eye burning	0.72*	0.72*
Overall eye dryness	0.71*	0.71*
Overall eye soreness	0.70*	0.70*
Overall eye stinging	0.69*	0.69*
Overall eye irritation	0.68*	0.68*
Overall eye redness	0.67*	0.67*
Overall eye itching	0.66*	0.66*
Overall eye watering	0.65*	0.65*
Overall eye burning	0.64*	0.64*
Overall eye dryness	0.63*	0.63*
Overall eye soreness	0.62*	0.62*
Overall eye stinging	0.61*	0.61*
Overall eye irritation	0.60*	0.60*
Overall eye redness	0.59*	0.59*
Overall eye itching	0.58*	0.58*
Overall eye watering	0.57*	0.57*
Overall eye burning	0.56*	0.56*
Overall eye dryness	0.55*	0.55*
Overall eye soreness	0.54*	0.54*
Overall eye stinging	0.53*	0.53*
Overall eye irritation	0.52*	0.52*
Overall eye redness	0.51*	0.51*
Overall eye itching	0.50*	0.50*
Overall eye watering	0.49*	0.49*
Overall eye burning	0.48*	0.48*
Overall eye dryness	0.47*	0.47*
Overall eye soreness	0.46*	0.46*
Overall eye stinging	0.45*	0.45*
Overall eye irritation	0.44*	0.44*
Overall eye redness	0.43*	0.43*
Overall eye itching	0.42*	0.42*
Overall eye watering	0.41*	0.41*
Overall eye burning	0.40*	0.40*
Overall eye dryness	0.39*	0.39*
Overall eye soreness	0.38*	0.38*
Overall eye stinging	0.37*	0.37*
Overall eye irritation	0.36*	0.36*
Overall eye redness	0.35*	0.35*
Overall eye itching	0.34*	0.34*
Overall eye watering	0.33*	0.33*
Overall eye burning	0.32*	0.32*
Overall eye dryness	0.31*	0.31*
Overall eye soreness	0.30*	0.30*
Overall eye stinging	0.29*	0.29*
Overall eye irritation	0.28*	0.28*
Overall eye redness	0.27*	0.27*
Overall eye itching	0.26*	0.26*
Overall eye watering	0.25*	0.25*
Overall eye burning	0.24*	0.24*
Overall eye dryness	0.23*	0.23*
Overall eye soreness	0.22*	0.22*
Overall eye stinging	0.21*	0.21*
Overall eye irritation	0.20*	0.20*
Overall eye redness	0.19*	0.19*
Overall eye itching	0.18*	0.18*
Overall eye watering	0.17*	0.17*
Overall eye burning	0.16*	0.16*
Overall eye dryness	0.15*	0.15*
Overall eye soreness	0.14*	0.14*
Overall eye stinging	0.13*	0.13*
Overall eye irritation	0.12*	0.12*
Overall eye redness	0.11*	0.11*
Overall eye itching	0.10*	0.10*
Overall eye watering	0.09*	0.09*
Overall eye burning	0.08*	0.08*
Overall eye dryness	0.07*	0.07*
Overall eye soreness	0.06*	0.06*
Overall eye stinging	0.05*	0.05*
Overall eye irritation	0.04*	0.04*
Overall eye redness	0.03*	0.03*
Overall eye itching	0.02*	0.02*
Overall eye watering	0.01*	0.01*
Overall eye burning	0.00*	0.00*

*Significant difference (p < 0.05) between Regimen 1 and comparison regimens.

RESULTS - Continued

Figure 2 - Subject Likert Questionnaire
"This solution gives my lenses long lasting comfort"




Figure 3 - Subject Likert Questionnaire
"This solution enhances my comfort"





Figure 4 - Visual Analog Scales



RESULTS - Continued

Figure 5 - Corneal Staining Severity




Figure 6 - Corneal Staining Severity





Figure 7 - Acute (1-day, 14 hours) & Long Term (30-day, 14 hours) Wetting Angle Measurements



RESULTS - Continued

Figure 8 - Mean Wetting Angle (29 Days, 14 hours)




Figure 9 - Correlation - Between Wetting Angle and Subjective Lens Moisture






Figure 10 - Mean Residual Lysozyme (Study 1 and Study 5)



RESULTS - Continued

Figure 11 - Rewetting Drop Use



Other Variables

- Average daily lens wear time ranged from 12.2 to 14.7 hours over the course of the five studies. No significant regimen differences were observed for average daily lens wear time.
- No clinically relevant changes from baseline in best-corrected Snellen visual acuity were observed.
- No clinically relevant differences in slit-lamp findings were observed between regimens during the studies.
- Adverse events related to any regimen were nonserious, mild or moderate in intensity, and resolved with or without treatment.

CONCLUSIONS

- Subjective comfort and regimen acceptability were consistently higher for Regimen 1 when compared with Regimens 2, 3, and 4, particularly in symptomatic subjects. Better comfort and acceptability could potentially lead to fewer contact lens dropouts.
- Corneal staining, both severity and area, was significantly greater for Regimen 2 and 3 when compared with Regimen 1. An uncompromised corneal epithelium is paramount to healthy lens wear.
- Contact lens wettability, assessed by ex vivo contact lens angle measurements, was significantly better for Regimen 1 compared with Regimens 2 and 3. Additionally, a moderate correlation between wettability and a subjective feeling of lens moisture was found.
- Residual lens lysozyme was significantly lower for Regimen 1 than for Regimens 2 and 3. Investigator-rated cleanliness was significantly better for subjects using Regimen 1 than for those using Regimen 2.
- Subjects using Regimen 2 or Regimen 4 used significantly more rewetting drops than those randomized to Regimen 1. The increased use of rewetting drops over the course of the study for some regimens might have diminished the ability to show comfort differences among the regimens in some instances.